


# Data Sheet

## Norit H2OK Medical TapFilter



### Dimensions & Weight

Overall length	136 mm
Overall diameter	Ø 88 mm
Dry weight system	240 grams

### Connection

To mains	22 mm female thread
Adapter	24 mm male thread

### Storage & Handling

- Keep dry during storage
- Protect against freezing after first use
- Handle with care; do not expose to shocks
- Disposable product; discard as regular waste

### Product name      Productcode

H2OK Medical TapFilter	31.2.201
Replacement Cartridge	31.2.251

### Packaging unit

Sterilized by gamma irradiation	
Starterset	12 pcs./box
Replacement cartridge	24 pcs./box



MEMBER



### Filtrix BV

Siliciumweg 28 • 3812 SX Amersfoort • The Netherlands

T +31 33 422 04 40 • F +31 33 422 04 41

E info@filtrix.com • I www.filtrix.com

Please visit our website to obtain information about your local support.

### Technology

Capillary microfiltration membranes

### Performance<sup>1</sup>

Rated capacity 2.000 liter  
Lifetime 31 days maximum  
Flowrate 6 l/min at 2 bar

### Compliance

Biological retention<sup>2</sup>      bacteria log 7 (B. Diminuta)  
Mechanical effects      turbidity reduction  
Material safety      extraction tests according to NSF/KTW

Structural integrity<sup>3</sup>      hydrostatic pressure 8 bar

### Operating conditions

Max. operating pressure      6 bar  
Operating temperature      min. 0 °C; max. 40 °C  
60 °C for 30 min. during lifetime  
Replacement interval      31 days (disposable product)  
Chlorine exposure      100 ppm

<sup>1</sup> Product performance is dependent on incoming water conditions

<sup>2</sup> Tested according to ASTM F838-05 by Vitens Laboratory, The Netherlands

<sup>3</sup> Tested according to NSF 42

**Limited warranty:** Filtrix BV represents and warrants for the warranty period of one year that its products will be free from substantial defects in materials and workmanship. This warranty is subject to exclusions and limitations. For details on warranty and limitation of liability please refer to our warranty sheet.

**Disclaimer:** The information and data contained in this document are based on our general experience and are believed to be correct. They are given in good faith and are intended to provide a guideline for the selection and use of our products. Since the conditions under which our products may be used are beyond our control, this information does not imply any guarantee of final product performance and we cannot accept any liability with respect to the use of our products. The quality of our products is guaranteed under our conditions of sale. Existing industrial property rights must be observed.

MFE E-10-2008

## Norit H2OK Medical TapFilter

**Norit**  
leading in purification

**Filtrix**